

Disposable Line


Craniotomy Cutter



The Craniotomy Cutter is suitable for use in conjunction with an electric or pneumatic craniotome to perform craniotomies and laminectomies.



Cutting Length




12 mm

16 mm

25 mm

Coupling Types



Technical Information



ANVISA REGISTER
80482490025
(ER-01 - ER-20)
80482499018
(ER-21 - ER-30)



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

