

Disposable Line


Round Fluted Bur



The Round Fluted Bur is recommended for aggressive grinding of bones and cartilage, with the purpose of removing osteocartilaginous tumors, performing decompressions in nervous or vascular structures, or prepare the area for the placement of bone fixators and implants.



Coupling Types



Technical Information



ANVISA REGISTER
80482490025
(ER-01 - ER-20)
80482499018
(ER-21 - ER-30)



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

