


## Disposable Line


### Round Diamond Bur



The Round Diamond Bur is indicated for delicate grinding and polishing of bone spicules, with the purpose of promoting bone finishing, it is also indicated for the decompression of nerve or vascular structures adjacent to the bone as it does not directly damage the soft tissues.



## Coupling Types



## Technical Information



ANVISA REGISTER  
80482490024  
(ER-01 - ER-20)  
80482499017  
(ER-21 - ER-30)



RAW MATERIAL  
STAINLESS STEEL



ESTERELIZAÇÃO  
ÓXIDO DE ETILENO



VALIDITY  
5 YEARS

