

Disposable Line

Lindemann Carbide Bur



The Lindemann Carbide Bur is suitable for performing osteotomies in bones and cartilage in places with tight access.

Carbide: The tungsten carbide stands out for its high hardness and resistance, providing precise and durable cuts, even in extremely hard materials. In the context of healthcare products, it offers several advantages, mainly in terms of cutting durability and thermal transfer.



diameter

1,6 mm




diameter

1,9 mm



Cutting Length



9,0 mm


15,2 mm

CUTTING LENGTH OPTIONS:

Ø1,6 mm: Cutting length 9,0 mm

Ø1,9 mm: Cutting length 15,2 mm

Coupling Types



ER-01



CARBIDE

Technical Information



ANVISA REGISTER
80482490028



RAW MATERIAL
CARBIDE



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS



www.biometal.med.br



+55 (19) 3495-2485



vendas@biometal.med.br