


Disposable Line

## Micro Sagittal Saw (BLS-01)



The Micro Sagittal Saw is often used in buccomaxillofacial surgery and small fragments to perform osteotomies with high precision and control through oscillations in the longitudinal direction.

**Model:**  
**BLS-01**



LENTH "A"	24,0 mm	27,0 mm	28,0 mm	30,0 mm	35,0 mm	40,0 mm	45,0 mm
WIDTH "B"	7,0 mm	9,0 mm					
THICKNESS	0,5 mm						

## Coupling Types



### Technical Information



ANVISA REGISTER  
80482499009



RAW MATERIAL  
STAINLESS STEEL



STERELIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS



[www.biometal.med.br](http://www.biometal.med.br)



+55 (19) 3495-2485



[vendas@biometal.med.br](mailto:vendas@biometal.med.br)