


Non-Disposable Line

Acetabular Flexible Drill




The Acetabular Flexible Drill is indicated for drilling the acetabular region in orthopedic hip replacement procedures (hip arthroplasty). The flexibility of the shaft allows surgeons to achieve specific drilling angles and positions, facilitating the successful insertion and fixation of implants.

CUTTING AREA LENGTH "A"	25 mm	30 mm	35 mm	40 mm	45 mm	50 mm	55 mm	60 mm
TOTAL LENGTH "B"	205 mm	210 mm	215 mm	220 mm	225 mm	230 mm	235 mm	240 mm



Ø diâmetro	Ø diâmetro	Ø diâmetro	Ø diâmetro
3,2 mm	3,5 mm	3,8 mm	4,5 mm

Coupling Types



Informações Técnicas



ANVISA REGISTER
80482499024



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
N/A



VALIDITY
N/A



www.biometal.med.br



(19) 3495-2485



vendas@biometal.med.br