


Disposable Line

Craniotomy Kit

! The Craniotomy Kit consists of a set of instruments designed to access the brain using the craniotomy procedure, one of the main objectives of which is to remove brain tumors, treat aneurysms or repair traumatic injuries to the brain.



Kit Contents

Ref	Qty	Descrição
01	01	Suction cannula
02	01	Cooling cannula
03	01	Craniotomy cutter
04	01	Twist drill (stop)
05	01	Auto stop cranial perforator
06	01	Silicone cooling hose (1,5m)

Technical Information



ANVISA REGISTER
80482490026



RAW MATERIAL
N/A



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

