

Disposable Line

# Round Fluted Bur



**!** The Round Fluted Bur is recommended for aggressive grinding of bones and cartilage, with the purpose of removing osteocartilaginous tumors, performing decompressions in nervous or vascular structures, or prepare the area for the placement of bone fixators and implants.

Ø diameter 0,6 mm	Ø diameter 0,8 mm	Ø diameter 1,0 mm	Ø diameter 1,4 mm	Ø diameter 2,3 mm	Ø diameter 3,1 mm	Ø diameter 4,0 mm
Ø diameter 4,5 mm	Ø diameter 5,0 mm	Ø diameter 6,0 mm	Ø diameter 7,0 mm	*Lenght upon request		



## Coupling Types

 ER-01	 ER-02	 ER-03	 ER-04	 ER-05
 ER-06	 ER-07	 ER-08	 ER-09	 ER-10
 ER-11	 ER-12	 ER-13	 ER-15	 ER-16
 ER-17	 ER-18	 ER-19	 ER-20	

### Technical Information



ANVISA REGISTER  
80482490025



RAW MATERIAL  
STAINLESS STEEL



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

