


Disposable Line

Round Diamond Bur



! The Round Diamond Bur is indicated for delicate grinding and polishing of bone spicules, with the purpose of promoting bone finishing, it is also indicated for the decompression of nerve or vascular structures adjacent to the bone as it does not directly damage the soft tissues.

0,6 mm	0,8 mm	1,0 mm	1,4 mm	2,3 mm	3,1 mm	4,0 mm	
4,5 mm	5,0 mm	6,0 mm	7,0 mm	*Lenght upon request			

Coupling Types

ER-01	ER-02	ER-03	ER-04	ER-05
ER-06	ER-07	ER-08	ER-09	ER-10
ER-11	ER-12	ER-13	ER-15	ER-16
ER-17	ER-18	ER-19	ER-20	

Technical Information



ANVISA REGISTER
80482490024



RAW MATERIAL
STAINLESS STEEL



ESTERELIZAÇÃO
ÓXIDO DE ETILENO



VALIDITY
5 YEARS

