


Disposable Line


Acorn Bur




The Acorn Bur is indicated for precise and safe bone grinding, its design provides less vibration and a smoother cut, it is generally used in neurosurgical, orthopedic and spine surgery.



*Lenght upon request



Coupling Types



Technical Information



ANVISA REGISTER
80482490025



RAW MATERIAL
STAINLESS STEEL



ESTERELIZAÇÃO
ÓXIDO DE ETILENO



VALIDITY
5 YEARS



www.biometal.med.br



+55 (19) 3495-2485



vendas@biometal.med.br