



Disposable Line


Match Head Diamond Bur

! The Match Head Diamond Bur is indicated for grinding and finishing of osteocartilaginous tissues in minimally invasive spinal procedures, foraminotomies, keyhole craniotomies, maxillocraniofacial and cranial base procedures.



EXTRA COARSE!

Coupling Types



Technical Information



ANVISA REGISTER
80482490024



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

