

Disposable Line

**Flexible Round Fluted SKB Bur (Minimal Invasive)**

**!** The Flexible Round Fluted SKB Bur is indicated for aggressive bone grinding and access to skull base structures such as tumor resections and hypophysis procedures.



∅ diameter  
**1,5 mm**

∅ diameter  
**2,0 mm**


∅ diameter  
**2,5 mm**








∅ diameter  
**3,0 mm**

∅ diameter  
**4,0 mm**

∅ diameter  
**5,0 mm**



**Coupling Types**



Technical Information



ANVISA REGISTER  
80482499005



RAW MATERIAL  
STAINLESS STEEL



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

