


Disposable Line

Match Head Diamond Skb Bur (Minimal Invasive)








! The flexible round diamond skb bur is indicated for aggressive abrasion and removal of bone spicules, in order to promote bone finishing, it can be used as an adjuvant in transfenoidal accesses where nerve and vascular structures are juxtaposed to the bone.



!
EXTRA COARSE!

Coupling Types



Technical Information



ANVISA REGISTER
80482499006



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

