

**Disposable Line**

**Bud Bur**




The Bud Bur is indicated for the controlled and progressive removal of bone and cartilage and is generally used for sinus ostia enlargement, mastoidectomies, osteoplasties and preparation of the area for implant placement.




∅ diameter  
**4,0 mm**

∅ diameter  
**5,0 mm**



Total Length




**45,0 mm**

**55,0 mm**

**70,0 mm**

**Coupling Types**



**Technical Information**



ANVISA REGISTER  
80482490025



RAW MATERIAL  
STAINLESS STEEL



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

