

Disposable Line

## Cross Cut Carbide Bur



The Cross Cut Carbide Bur is indicated for making perforations, precise cuts and delicate osteotomies in bones and cartilage, as well as being effective in dental removals.



699



701






702



703



Cutting Length




**CUTTING LENGTH OPTIONS:**

699: Cutting length 4,6mm

701: Cutting length 5,1mm

702 / 703: Cutting length 5,1mm or 7,6mm

## Coupling Types



### Technical Information



ANVISA REGISTER  
80482490028



RAW MATERIAL  
CARBIDE



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

