




Disposable Line

## Lindemann Carbide Bur



The Lindemann Carbide Bur is suitable for performing osteotomies in bones and cartilage in places with tight access.




**CUTTING LENGTH OPTIONS:**

∅1,6 mm: Cutting length 9,0 mm

∅1,9 mm: Cutting length 15,2 mm

## Coupling Types



### Technical Information



ANVISA REGISTER  
80482490028



RAW MATERIAL  
CARBIDE



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

