

Disposable Line


# Micro Sagittal Saw (BLS-01)



The Micro Sagittal Saw is often used in buccomaxillofacial surgery and small fragments to perform osteotomies with high precision and control through oscillations in the longitudinal direction.



Model:  
**BLS-01**



LENTH "A"	24,0 mm	27,0 mm	28,0 mm	30,0 mm	35,0 mm	40,0 mm	45,0 mm
WIDTH "B"	7,0 mm		9,0 mm				
THICKNESS	0,5 mm						

## Coupling Types



### Technical Information



ANVISA REGISTER  
80482499009



RAW MATERIAL  
STAINLESS STEEL



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

