



Disposable Line

# Micro Reciprocating Saw (BLR-01)




**!** The Micro Reciprocating Saw is often used in oral and maxillofacial surgery and small fragments to perform osteotomies with high precision and control through oscillations in the vertical direction.

**Model:**  
**BLR-01**



<b>CUTTING AREA "A"</b>	22 mm					27 mm				
<b>SHAFT LENGHT "B"</b>	29,5 mm	40,0mm	42,7 mm	56,0 mm	60,0 mm	29,5 mm	40,0mm	42,7 mm	56,0 mm	60,0 mm
<b>SAW THICKNESS</b>	0,4 mm									

## Coupling Types



### Technical Information



ANVISA REGISTER  
80482499009



RAW MATERIAL  
STAINLESS STEEL



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

