


Disposable Line


Micro Reciprocating Saw (BLR-02)



The Micro Reciprocating Saw is often used in oral and maxillofacial surgery and small fragments to perform osteotomies with high precision and control through oscillations in the vertical direction.



Modelo:
BLR-02



CUTTING AREA "A"	10,0 mm	20,0mm	30,0 mm
SAW THICKNESS	0,3 mm		

Coupling Types



Technical Information



ANVISA REGISTER
80482499009



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

