

Disposable Line

Percutaneous Burs (Hand and Foot)

! The Percutaneous Burs are indicated for minimal invasive procedures on small fragments with the purpose of drilling, cutting, grinding, finishing and bone scraping.




Model

Description and Dimension

	Straight linear cutter for percutaneous osteotomy ø2,0mm x 13mm
	Straight linear cutter for percutaneous osteotomy ø2,2mm x 25mm
	Straight helical cutter for percutaneous osteotomy ø2,0mm x 8mm
	Straight helical cutter for percutaneous osteotomy ø2,2mm x 12mm
	Straight helical cutter for percutaneous osteotomy ø2,9mm x 13,2mm
	Conical helical cutter for percutaneous osteotomy ø4,3mm x 13mm
	Conical helical cutter for percutaneous osteotomy ø5,0mm x 15mm.

Coupling Types



Technical Information



ANVISA REGISTER
80482490025



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

