



Non-Disposable Line

Orthopedic Flexible Reamer




The Orthopedic Flexible Reamer is indicated for the creation of intramedullary canals in long bones, such as the femur and tibia. This tool has a flexible shaft to facilitate access to anatomically complex areas and allows the removal of bone efficiently, with the aim of facilitating the insertion of medical devices such as intramedullary rods to stabilize fractures and promote proper bone regeneration.

- 250 mm
- 275 mm
- 300 mm
- 325 mm
- 350 mm
- 375 mm
- 400 mm
- 450 mm
- 475 mm
- 500 mm



- | | | | | | | |
|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| ∅ diameter
6,0 mm | ∅ diameter
6,5 mm | ∅ diameter
7,0 mm | ∅ diameter
7,5 mm | ∅ diameter
8,0 mm | ∅ diameter
8,5 mm | ∅ diameter
9,0 mm |
| ∅ diameter
9,5 mm | ∅ diameter
10,0 mm | ∅ diameter
10,5 mm | ∅ diameter
11,0 mm | ∅ diameter
11,5 mm | ∅ diameter
12,0 mm | ∅ diameter
12,5 mm |
| ∅ diameter
13,0 mm | ∅ diameter
13,5 mm | ∅ diameter
14,0 mm | ∅ diameter
14,5 mm | ∅ diameter
15,0 mm | ∅ diameter
15,5 mm | ∅ diameter
16,0 mm |

Coupling Types



Technical Information



ANVISA REGISTER
80482490019



RAW MATERIAL
STAINLESS STEEL



ESTERILIZAÇÃO
N/A



VALIDITY
N/A

