

Non-Disposable Line

Acetabular Flexible Drill




The Acetabular Flexible Drill is indicated for drilling the acetabular region in orthopedic hip replacement procedures (hip arthroplasty). The flexibility of the shaft allows surgeons to achieve specific drilling angles and positions, facilitating the successful insertion and fixation of implants.




- 250 mm
- 260 mm
- 270 mm
- 280 mm
- 290 mm
- 300 mm
- 310 mm
- 320 mm

Total Length



- ϕ
diameter
3,2 mm
- ϕ
diameter
4,5 mm

Coupling Types



Technical Information



ANVISA REGISTER
80482490019



RAW MATERIAL
STAINLESS STEEL



ESTERELIZAÇÃO
N/A



VALIDITY
N/A

