


Disposable Line

Micro Sagital Saw



The Micro Sagital Saw is often used in buccomaxillofacial surgery and small fragments to perform osteotomies with high precision and control through oscillations in the longitudinal direction.

Model:
BLS-02



LENTH "A"	25,0 mm		
WIDTH "B"	6,0 mm	8,0 mm	10,0 mm
THICKNESS	0,3 mm		

Coupling Types



Technical Information



ANVISA REGISTER
80482499009



RAW MATERIAL
STAINLESS STEEL



STERELIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS



www.biometal.med.br



(19) 3495-2485



vendas@biometal.med.br