


Disposable Line

Micro Basal Reciprocating Saw (Right / Left)



The Basal Reciprocating Saw is indicated in bucomaxillofacial surgeries for procedures requiring controlled and precise bone cuts, minimizing damage to adjacent structures. This blade is particularly useful as it allows access to complex anatomical areas.



SAW WIDTH
"A"

9,85 mm



ROD LENTH
"B"

32,60 mm

ROD LENTH
"C"

44,20 mm

Coupling Types



Technical Information



ANVISA REGISTER
80482499009



RAW MATERIAL
STAINLESS STEEL



STERELIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS



www.biometal.med.br



(19) 3495-2485



vendas@biometal.med.br