



Disposable Line

Micro Oscilating Saw





! The Micro Oscilating Saw is often used in buccomaxillofacial surgery to perform osteotomies with high precision and control through oscillations in the longitudinal direction, especially in areas that are difficult to access.



SAW WIDTH "A"	SAW LENTH "B"	THICKNESS "C"	TOTAL LENGHT "C"
6,0 mm	10 mm 15 mm	0,4 mm	30 mm 45 mm
6,5 mm			
9,2 mm			

Coupling Types



Technical Information



ANVISA REGISTER
80482499009



RAW MATERIAL
STAINLESS STEEL



STERELIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS



www.biometal.med.br



(19) 3495-2485



vendas@biometal.med.br