



Disposable Line

Cutting Debridement Cannula




The cutting debridement cannula is used in minimally invasive surgeries, particularly in spinal arthroscopy, for removing, cutting, drilling, and scraping cartilage and bone tissue, facilitating decompression, debridement, shaping, and precision in the surgical field.

Models:



A	Outer Diameter	3,5 mm		
B	Length	300 mm	330 mm	360 mm

Coupling Types



Technical Information



ANVISA REGISTER
80482499025



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS

