


Disposable Line

# Diamond Debridement Cannula



The diamond debridement cannula is used in minimally invasive surgeries, particularly in spinal arthroscopy, for scraping, polishing, and finishing cartilage and bone tissue, facilitating decompression, debridement, shaping, and precision in the surgical field.


## Models:




**BCDD  
001**



**BCDD  
002**



**BCDDPL  
001**




**BCDDPL  
002**



**BCDDPF  
001**





**BCDDPF  
002**



<b>A</b>	Outer Diameter	3,5 mm		
<b>B</b>	Length	300 mm	330 mm	360 mm

## Coupling Types



## Technical Information



ANVISA REGISTER  
80482499026



RAW MATERIAL  
STAINLESS STEEL



STERILIZATION  
ETHYLENE OXIDE



VALIDITY  
5 YEARS

