



Disposable Line

Cutting Microdebridement Cannula





The cutting microdebridement cannula is used in minimally invasive surgeries, mainly in the areas of orthopedics and otolaryngology, to remove, cut, and/or debride soft or cartilaginous tissue in a precise and controlled manner.

Models:



A	Outer Diameter	4,5 mm
B	Length	95 mm 115 mm

Coupling Types



Technical Information



ANVISA REGISTER
80482499031



RAW MATERIAL
STAINLESS STEEL



STERILIZATION
ETHYLENE OXIDE



VALIDITY
5 YEARS



www.biometal.med.br



(19) 3495-2485



vendas@biometal.med.br